


EQUIPARABILI

Tray & Implant instruments




English version




MADE IN ITALY






Implant prosthetic surgery, in recent years, has made possible fixed prosthetic rehabilitation solutions that have become essential in daily dental practice.

Even a dentist today, to meet the modern needs of the patient, must take professional action and equip themselves with appropriate technical tools to achieve these goals.

Jetimplant Srl, with its kits, aims to provide dentists with modern tools that allow them to save time, improve work ergonomics, and simplify surgery.

The fundamental characteristics of the comparable implantology and surgery kits can be described as follows:

- The anodized aluminum kits consist of a base and a lid; inside, there are two inserts. One is specifically dedicated to the surgical procedure, where the measurements of the various instruments are described; the second insert is equipped with retention rubber pads and functions to collect all the surgical instruments selected by the operator.
- The base, the lid and the two inserts are assembled together with interlocking. This allows the operator easy cleaning, disinfection, sterilization, and reassembly.
- In the comparable kits, the instrument pickup point is devoid of rubber pads or other retention elements on the surface, thus facilitating cleaning and washing after being contaminated by the operator. The underlying insert, where the retention rubber pads are housed, is therefore protected by the upper insert and kept away from contact.
- The comparable kits, due to their particular design, allow for simple packing and sterilization maneuvers in a steam autoclave.
- All surgical instruments, once placed in their designated slots, cannot be moved or overturned in any way during handling.
- The comparable kits are opened in the operating field, facilitating the surgeon in choosing and using the instruments (drills and surgical accessories).
- **The comparable kits are protected by two types of patents.**



Kit Identikit



KIT Identikit 1

Contents

DESCRIPTION

- 11 Identikit threaded instruments
- 2 Identikit threaded instruments IMZ®
- 13 Identikit hexagonal instruments
- 1 Identikit square instrument
- 1 Identikit octagonal instrument

Kit Code: JETID1


KIT Identikit 2

Contents

DESCRIPTION

- 1 Torx Key 5,5
- 1 Torx Key 6,0
- 1 Hexagonal Key 0,9
- 1 Hexagonal Key 1,15
- 1 Hexagonal Key 1,2
- 1 Hexagonal Key 1,25
- 1 Hexagonal Key 1,3
- 1 Hexagonal Key 1,5
- 1 Hexagonal Key 1,7
- 1 Square Key 0,9
- 1 Slot Key 1,6
- 1 Locator Key
- 1 Equator Key
- 1 TPA Key
- 1 Multy Unit Key
- 1 Ratchet
- 11 Identikit threaded instruments
- 2 Identikit threaded instruments IMZ®
- 13 Identikit hexagonal instruments
- 1 Identikit square instrument
- 1 Identikit octagonal instrument

Kit Code: JETID2




Kit 118

Standard

Contents

The kit consists of 15 screwdrivers:

- | | |
|---------------|--------------|
| 01 · Torx 5,5 | 01 · Torx 6 |
| 01 · E 0,9 | 01 · E 1,15 |
| 01 · E 1,2 | 01 · E 1,25 |
| 01 · E 1,3 | 01 · E 1,5 |
| 01 · E 1,7 | 01 · Q 0,9 |
| 01 · T 1,6 | 01 · LOCATOR |
| 01 · EQUATOR | 01 · TPA 6 |
| 01 · MUA | |




Kit 118

Medium

Contents

The kit consists of 15 screwdrivers:

01 · Torx 5,5	01 · Torx 6
01 · E 0,9	01 · E 1,15
01 · E 1,2	01 · E 1,25
01 · E 1,3	01 · E 1,5
01 · E 1,7	01 · Q 0,9
01 · T 1,6	01 · LOCATOR
01 · EQUATOR	01 · TPA 6
01 · MUA	




Kit 118 Micromotor

Medium

Contents

The kit consists of 16 screwdrivers:

01 · Torx 4	01 · Torx 5,5
01 · Torx 6	01 · E 0,9
01 · E 1,15	01 · E 1,2
01 · E 1,25	01 · E 1,27
01 · E 1,3	01 · E 1,5
01 · E 1,7	01 · Slot 1,6
01 · TPA 4	01 · TPA 5
01 · TPA 6	01 · TPA STR




Kit 118 Micromotor

Small

Contents

The kit consists of 12 screwdrivers:

01 · Torx 4	01 · Tork 5,5
01 · Tork 6	01 · E 0,9
01 · E 1	01 · E 1,15
01 · E 1,2	01 · E 1,25
01 · E 1,3	01 · E 1,5
01 · E 1,7	01 · Slot 1,6




Kit TPA for Angled Hole Screws Micromotor

Medium

Contents

The kit consists of:

- 1 manual handle
- 1 short manual knob
- 1 TPA 4-lobed screwdriver
- 1 TPA 5-lobed screwdriver
- 1 TPA 6-lobed screwdriver
- 1 Straumann®-type TPA




Kit TPA for Angled Hole Screws Micromotor

Small

Contents

The kit consists of:

- 1 · 6-lobed screwdriver
- 1 · 5-lobed screwdriver
- 1 · 4-lobed screwdriver
- 1 · STR hexagonal screwdriver




Kit Mucotomes

Small

Contents

The kit consists of 5 mucotomes:

- Mucotome Ø 3mm
- Mucotome Ø 4mm
- Mucotome Ø 4.5mm
- Mucotome Ø 5mm
- Mucotome Ø 5.5mm




Kit Extractors

Medium

Contents

The kit consists of:

- 1 ratchet
(torque wrench)
- 4 extractors
(2 for through screws
with worn hexagon;
2 for fractured implant)



Kit Laboratory Screwdriver

Standard



Contents


The kit consists of 6 screwdrivers:

Torx 5,5
Torx 6
E 1,2
E 1,25
1 EQUATOR
1 LOCATOR

Kit Laboratory Screwdriver

TPA

Medium



Contents

The kit consists of 4 screwdrivers:

TPA 4 LOBI
TPA 5 LOBI
TPA 6 LOBI
TPA STR

Kit Laboratory Screwdriver

Medium



Contents


The kit consists of 6 screwdrivers:

Torx 5,5
Torx 6
E 1,2
E 1,25
1 EQUATOR
1 LOCATOR

Kit Laboratory Screwdriver

TPA

Small



Contents

The kit consists of 4 screwdrivers:

TPA 4 LOBI
TPA 5 LOBI
TPA 6 LOBI
TPA STR

Kit Laboratory Screwdriver

Small



Contents

The kit consists of 6 screwdrivers:

Torx 5,5
Torx 6
E 1,2
E 1,25
1 EQUATOR
1 LOCATOR



JETIMPLANT S.r.l.

Via privata Emilia Romagna, 13/15
26855 Lodi Vecchio (LO), Italy
Telephone: (+39) 0371 46.03.02
Whatsapp: (+39) 351 46.14.739
www.jetimplant.com
jetimplantsrl@gmail.com